



Here comes the sun. Upgrade from XHL to the LED^{HQ} module.

CHANGING BULBS IS NOW A THING OF THE PAST!

With this new option, the sun is rising for many owners of HEINE instruments with XHL illumination.

- ▶ Virtually unlimited durability with LED^{HQ} in HEINE quality
- ▶ Maintenance-free and reliably ready for use
- ▶ Longer running times for both single-use batteries and rechargeable batteries thanks to lower power consumption
- ▶ Sustainable, as this is not a single-use product and there is less battery waste
- ▶ Tested spare part guarantees patient safety*

| Instrument | Type | Power source | Previous XHL bulb | New LED ^{HQ} module | New order code |
|----------------|--|--------------------------|-------------------|------------------------------|----------------|
| Otoscope | BETA 400 F.O. BETA 200 F.O. K180 | Battery handle 2,5V | X-001.88.077 | X-028.88.132 | 1610326 |
| | | Rechargeable handle 3,5V | X-002.88.078 | X-008.88.133 | 1610327 |
| Ophthalmoscope | BETA 200S BETA 200 | Battery handle 2,5V | X-001.88.069 | X-028.88.134 | 1610328 |
| | | Rechargeable handle 3,5V | X-002.88.070 | X-008.88.135 | 1610329 |
| Retinoscope | BETA 200 | Battery handle 2,5V | X-001.88.087 | X-028.88.130 | 1610421 |
| | | Rechargeable handle 3,5V | X-002.88.089 | X-008.88.131 | 1610330 |
| Veterinary | G100 slit illumination head BETA 200 VET otoscope | Battery handle 2,5V | X-001.88.077 | X-028.88.132 | 1610326 |
| | | Rechargeable handle 3,5V | X-002.88.078 | X-008.88.133 | 1610327 |

*Only the original is as good as the original. Use HEINE ORIGINAL PARTS and not third-party products from other manufacturers. Only then will our guarantee apply and performance features such as colour rendering, function and safety of our examination instruments are always guaranteed.