



APPLANATION TONOMETER AT 900 & AT 900 D

Go for the Gold Standard

Tradition and Innovation – Since 1858, visionary thinking and a fascination with technology have guided us to develop innovative products of outstanding reliability: Anticipating trends to improve the quality of life.

HS HAAG-STREIT
DIAGNOSTICS

AT 900 and AT 900 D*

The gold standard in IOP measurement

Because it is fast and reproducible, Goldmann applanation tonometry has been the gold standard in tonometry for many decades. It is based on the Goldmann principle – tried and tested during more than 50 years of usage with millions of patients.

Furthermore, thanks to the high precision of manufacturing at Haag-Streit, the Goldmann applanation tonometer offers reliable performance over many years without the need for manufacturer maintenance.

AT 900 is the original version of the Goldmann tonometer and works 100% mechanically. The digital Goldmann tonometer AT 900 D* adds the advantages of state-of-the-art digital technology.

*AT 900 D is not available in the US and Canada



Swiss made

TONOSAFE – THE DISPOSABLE TONOMETER PRISMS

Convenient – Safe – Accurate

Tonosafe avoids the need for time-consuming and laborious cleaning and disinfection without compromising the original Goldmann quality. It has been developed for single use only and thus reliably eliminates the risk of cross-infection even in the most difficult situations. Calibrated to Goldmann standards, Tonosafe delivers accurate results with all Goldmann and Perkins contact tonometers.



TONOSAFE






GOLDMANN TONOMETER PRISMS

GOLDMANN TONOMETER PRISMS

Excellent optical quality

Goldmann measuring prisms are distinguished by their optical and mechanical performance which allows for fast and reliable measurements in daily practice. They are largely handmade by experienced workers, and are checked for tight tolerances in weight, optical quality, durability and for consistent measurement results.

Slit lamp model	mod. R	mod. T	mod. BQ
BD 900	■		
BM 900	■	■	
BI 900	■	■	
BP 900		■	■
BQ 900		■	■



Members of HAAG-STREIT Group

HAAG-STREIT Holding AG
www.haag-streit-holding.com

CLEMENT CLARKE Ltd.
www.clement-clarke.com

HAAG-STREIT AG, Diagnostics
www.haag-streit.com

HS DOMS GmbH
www.hs-doms.com

HAAG-STREIT AG, Verkauf Schweiz
www.haag-streit.ch

I PRO GmbH
www.ipro.de

HAAG-STREIT Deutschland GmbH
www.haag-streit.de

John Weiss Ltd.
www.johnweiss.com

HAAG-STREIT Far East
www.haag-streit-fareast.com

Möller-Wedel GmbH & Co KG
www.haag-streit-surgical.com

HAAG-STREIT Medtech AG
www.haag-streit-medtech.com

Möller-Wedel Optical GmbH
www.moeller-wedel-optical.com

HAAG-STREIT Surgical GmbH
www.haag-streit-surgical.com

OptoMedical Technologies GmbH
www.haag-streit-surgical.com

HAAG-STREIT UK
www.haag-streit-uk.com

Reliance Medical Inc.
www.haag-streit-usa.com

HAAG-STREIT USA
www.haag-streit-usa.com

SPECTROS AG
www.spectros.ch

Distributed by

Van Hopplynus Ophtalm

Rijksweg 10 - 2880 Bornem
Call toll free (Belgium) 0800 76 773
www.vho.be - info@vho.be